

For Research Use Only

CSN1S1 Polyclonal antibody

Catalog Number: 22648-1-AP

1 Publications



Basic Information

Catalog Number:

22648-1-AP

Size:

900 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18529

GenBank Accession Number:

BC128227

GeneID (NCBI):

1446

UNIPROT ID:

P47710

Full Name:

casein alpha s1

Calculated MW:

185 aa, 22 kDa

Observed MW:

21-24 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : human milk,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Suhua Chen	35565755	Nutrients	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human milk were subjected to SDS PAGE followed by western blot with 22648-1-AP (CSN1S1 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours.