

For Research Use Only

SRC-3/AIB1 Polyclonal antibody

Catalog Number: 22707-1-AP



Basic Information

Catalog Number:

22707-1-AP

Size:

280 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18424

GenBank Accession Number:

BC119001

GeneID (NCBI):

8202

UNIPROT ID:

Q9Y6Q9

Full Name:

nuclear receptor coactivator 3

Calculated MW:

1424 aa, 155 kDa

Observed MW:

155 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

FC, WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Raji cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Raji cells were subjected to SDS PAGE followed by western blot with 22707-1-AP (SRC-3/AIB1 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.



1×10^6 Raji cells were stained with 0.2ug SRC-3/AIB1 antibody (22707-1-AP, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.