

For Research Use Only

CHD5 Polyclonal antibody

Catalog Number: 22725-1-AP



Basic Information

Catalog Number:

22725-1-AP

Concentration:

800 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18621

GenBank Accession Number:

BC148803

GeneID (NCBI):

26038

UNIPROT ID:

Q8TD10

Full Name:

chromodomain helicase DNA binding protein 5

Calculated MW:

1954 aa, 223 kDa

Observed MW:

150-160 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse brain tissue, rat brain tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 22725-1-AP (CHD5 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.