

For Research Use Only

CTAG1B Polyclonal antibody

Catalog Number: 22785-1-AP



Basic Information

Catalog Number:

22785-1-AP

Size:

700 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18784

GenBank Accession Number:

BC130362

GeneID (NCBI):

1485

UNIPROT ID:

P78358

Full Name:

cancer/testis antigen 1B

Calculated MW:

180 aa, 18 kDa

Observed MW:

20 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HT-1080 cells, human testis tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 22785-1-AP (CTAG1B antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.