

For Research Use Only

# CTAG1B Polyclonal antibody, PBS Only

Catalog Number:22785-1-PBS



## Basic Information

Catalog Number:

22785-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18784

GenBank Accession Number:

BC130362

GeneID (NCBI):

1485

UNIPROT ID:

P78358

Full Name:

cancer/testis antigen 1B

Calculated MW:

180 aa, 18 kDa

Observed MW:

20 kDa

Purification Method:

Antigen Affinity purified

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 22785-1-AP (CTAG1B antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 22785-1-PBS in a different storage buffer formulation.