

For Research Use Only

# PYGO1 Polyclonal antibody

Catalog Number: 22799-1-AP

2 Publications



## Basic Information

Catalog Number:

22799-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18814

GenBank Accession Number:

BC146572

GeneID (NCBI):

26108

UNIPROT ID:

Q9Y3Y4

Full Name:

pygopus homolog 1 (Drosophila)

Calculated MW:

419 aa, 45 kDa

Observed MW:

45 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human, mouse

Positive Controls:

WB: HeLa cells,

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Jie Shen	36085540	Hum Cell	WB
Li Lin	33635162	Am J Physiol Heart Circ Physiol	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 22799-1-AP (PYGO1 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.