

For Research Use Only

PYGO1 Polyclonal antibody

Catalog Number: 22799-1-AP **2 Publications**



Basic Information

Catalog Number: 22799-1-AP	GenBank Accession Number: BC146572	Purification Method: Antigen affinity purification
Size: 300 µg/ml	GeneID (NCBI): 26108	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q9Y3Y4	
Isotype: IgG	Full Name: pygopus homolog 1 (Drosophila)	
Immunogen Catalog Number: AG18814	Calculated MW: 419 aa, 45 kDa	
	Observed MW: 45 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells,
Cited Applications: WB	
Species Specificity: human	
Cited Species: human, mouse	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Jie Shen	36085540	Hum Cell	WB
Li Lin	33635162	Am J Physiol Heart Circ Physiol	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 22799-1-AP (PYGO 1 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.