

For Research Use Only

Cytokeratin 28 Polyclonal antibody

Catalog Number: 22813-1-AP



Basic Information

Catalog Number:

22813-1-AP

Concentration:

550 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18932

GenBank Accession Number:

BC148794

GeneID (NCBI):

162605

UNIPROT ID:

Q7Z3Y7

Full Name:

keratin 28

Calculated MW:

464 aa, 51 kDa

Observed MW:

51-53 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse skin tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 22813-1-AP (KRT28 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.