

For Research Use Only

# BMP10 Polyclonal antibody

Catalog Number: 22858-1-AP



## Basic Information

Catalog Number:

22858-1-AP

Size:

150 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18746

GenBank Accession Number:

BC105063

GeneID (NCBI):

27302

UNIPROT ID:

O95393

Full Name:

bone morphogenetic protein 10

Calculated MW:

424 aa, 48 kDa

Observed MW:

48 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells, MCF7 cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



MCF7 cells were subjected to SDS PAGE followed by western blot with 22858-1-AP (BMP10 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.



HEK-293 cells were subjected to SDS PAGE followed by western blot with 22858-1-AP (BMP10 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.