

For Research Use Only

CDX4 Polyclonal antibody

Catalog Number: 22943-1-AP



Basic Information

Catalog Number:

22943-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19078

GenBank Accession Number:

BC128233

GeneID (NCBI):

1046

UNIPROT ID:

O14627

Full Name:

caudal type homeobox 4

Calculated MW:

284 aa, 30 kDa

Observed MW:

37 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : SKOV-3 cells, mouse thymus tissue, Sp2/O cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



SKOV-3 cells were subjected to SDS PAGE followed by western blot with 22943-1-AP (CDX4 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.