

For Research Use Only

BARD1 Polyclonal antibody

Catalog Number: 22964-1-AP

3 Publications



Basic Information

Catalog Number:

22964-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18965

GenBank Accession Number:

BC126426

GeneID (NCBI):

580

UNIPROT ID:

Q99728

Full Name:

BRCA1 associated RING domain 1

Calculated MW:

777 aa, 87 kDa

Observed MW:

90-100 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : SH-SY5Y cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yixi Xu	29106372	Elife	WB
Gozde Inci	36500714	Molecules	WB
Yunshang Chen	38296976	Cell Death Dis	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with 22964-1-AP (BARD1 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.