

For Research Use Only

DSG4 Polyclonal antibody

Catalog Number: 23087-1-AP



Basic Information

Catalog Number: 23087-1-AP	GenBank Accession Number: BC132907	Purification Method: Antigen Affinity purified
Size: 1000 µg/ml	GeneID (NCBI): 147409	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q865J6	
Isotype: IgG	Full Name: desmoglein 4	
Immunogen Catalog Number: AG19520	Calculated MW: 1040 aa, 114 kDa	
	Observed MW: 160 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse skin tissue,
Species Specificity: human, mouse	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse skin tissue were subjected to SDS PAGE followed by western blot with 23087-1-AP (DSG4 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.