

For Research Use Only

SCT Polyclonal antibody

Catalog Number: 23125-1-AP **1 Publications**



Basic Information

Catalog Number: 23125-1-AP	GenBank Accession Number: BC146571	Purification Method: Antigen Affinity purified
Size: 550 µg/ml	GeneID (NCBI): 6343	Recommended Dilutions: IHC 1:20-1:200
Source: Rabbit	UNIPROT ID: P09683	
Isotype: IgG	Full Name: SCT	
Immunogen Catalog Number: AG19419	Calculated MW: 121 aa, 13 kDa	

Applications

Tested Applications: IHC, ELISA	Positive Controls: IHC : human liver cancer tissue, human stomach tissue
Cited Applications: IF	
Species Specificity: human	
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Hirohumi Suzuki	25458814	Tissue Cell	IF

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human liver cancer tissue slide using 23125-1-AP (Secretin antibody at dilution of 1:50 (under 40x lens).



Immunohistochemical analysis of paraffin-embedded human stomach tissue slide using 23125-1-AP (Secretin antibody at dilution of 1:50 (under 40x lens).