

For Research Use Only

# CADM4 Polyclonal antibody, PBS Only

Catalog Number: 23136-1-PBS



## Basic Information

Catalog Number:

23136-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18922

GenBank Accession Number:

BC141627

GeneID (NCBI):

199731

UNIPROT ID:

Q8NFZ8

Full Name:

cell adhesion molecule 4

Calculated MW:

388 aa, 43 kDa

Observed MW:

65-70 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

## Background Information

CADM4, also known as IGSF4C or TSL2, is a cell adhesion molecule belonging to the immunoglobulin superfamily. CADM4 has been shown to serve as a tumor suppressor.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



RAW264.7 cells were subjected to SDS PAGE followed by western blot with 23136-1-AP (CADM4 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 23136-1-PBS in a different storage buffer formulation.