

For Research Use Only

ZFR Polyclonal antibody

Catalog Number: 23174-1-AP



Basic Information

Catalog Number: 23174-1-AP	GenBank Accession Number: BC137085	Purification Method: Antigen Affinity purified
Size: 400 µg/ml	GeneID (NCBI): 51663	Recommended Dilutions: WB 1:500-1:2000 IHC 1:50-1:500
Source: Rabbit	UNIPROT ID: Q96KR1	
Isotype: IgG	Full Name: zinc finger RNA binding protein	
Immunogen Catalog Number: AG19559	Calculated MW: 1074 aa, 117 kDa	
	Observed MW: 120 kDa	

Applications

Tested Applications:

IHC, WB, ELISA

Species Specificity:

human

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : HEK-293 cells,

IHC : human small intestine tissue, human heart tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

