

For Research Use Only

DLX6 Polyclonal antibody, PBS Only

Catalog Number: 23216-1-PBS

Featured Product



Basic Information

Catalog Number:

23216-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19724

GenBank Accession Number:

BC069363

GeneID (NCBI):

1750

UNIPROT ID:

P56179

Full Name:

distal-less homeobox 6

Calculated MW:

175 aa, 20 kDa

Observed MW:

32 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



NCCIT cell were subjected to SDS PAGE followed by western blot with 23216-1-AP (DLX6 Antibody) at dilution of 1:100 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 23216-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse embryo tissue slide using 23216-1-AP (DLX6 Antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 23216-1-PBS in a different storage buffer formulation.