

For Research Use Only

TBX18 Polyclonal antibody

Catalog Number: 23237-1-AP **2 Publications**



Basic Information

Catalog Number: 23237-1-AP	GenBank Accession Number: BC157841	Purification Method: Antigen Affinity purified
Size: 600 µg/ml	GeneID (NCBI): 9096	Recommended Dilutions: WB 1:500-1:2000 IHC 1:20-1:200
Source: Rabbit	UNIPROT ID: O95935	
Isotype: IgG	Full Name: T-box 18	
Immunogen Catalog Number: AG19258	Calculated MW: 607 aa, 65 kDa	
	Observed MW: 65-70 kDa	

Applications

Tested Applications:
IHC, WB, ELISA

Cited Applications:
WB, IF

Species Specificity:
human, mouse, rat

Cited Species:
human, rat

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB: HeLa cells, HEK-293 cells, human placenta tissue, mouse kidney tissue, rat kidney tissue

IHC: human kidney tissue, human breast cancer tissue

Background Information

The T-box (Tbx) motif is present in a family of genes whose structural features and expression patterns support their involvement in developmental gene regulation. The Tbx gene family are largely conserved throughout metazoan evolution, and these genes code for putative transcription factors that share a uniquely defining DNA-binding domain. Tbx genes are a family of developmental regulators with more than 20 members recently identified in invertebrates and vertebrates. Mutations in Tbx genes are associated with the onset of several human diseases. Our understanding of functional mechanisms of Tbx products has come mainly from the prototypical T/Brachyury, which is a transcription activator. The Tbx genes constitute a family of transcriptional regulatory genes that are implicated in a variety of developmental processes ranging from the formation of germ layers to the organizational patterning of the central nervous system. This antibody specifically reacts with the 65kd TBX18 protein.

Notable Publications

Author	Pubmed ID	Journal	Application
Xiuya Li	34488017	EBioMedicine	IF
Pei Xu	36484369	Mol Med Rep	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 23237-1-AP (TBX18 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human breast cancer tissue slide using 23237-1-AP (TBX18 Antibody) at dilution of 1:50 (under 40x lens).



Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using 23237-1-AP (TBX18 Antibody) at dilution of 1:50 (under 40x lens).