

For Research Use Only

FBXL19 Polyclonal antibody

Catalog Number: 23624-1-AP



Basic Information

Catalog Number:

23624-1-AP

Size:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20198

GenBank Accession Number:

BC059173

GeneID (NCBI):

54620

UNIPROT ID:

Q6PCT2

Full Name:

F-box and leucine-rich repeat protein

19

Calculated MW:

694 aa, 76 kDa

Observed MW:

60-75 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells, HeLa cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 23624-1-AP (FBXL19 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.