

For Research Use Only

ENOX1 Polyclonal antibody

Catalog Number: 23750-1-AP



Basic Information

Catalog Number:

23750-1-AP

Size:

900 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20715

GenBank Accession Number:

BC024178

GeneID (NCBI):

55068

UNIPROT ID:

Q8TC92

Full Name:

ecto-NOX disulfide-thiol exchanger 1

Calculated MW:

73 kDa

Observed MW:

73 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Jurkat cells, MCF-7 cells, Raji cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 23750-1-AP (ENOX1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.