

For Research Use Only

ZNF462 Polyclonal antibody

Catalog Number: 23793-1-AP



Basic Information

Catalog Number:

23793-1-AP

Concentration:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20699

GenBank Accession Number:

BC036884

GeneID (NCBI):

58499

UNIPROT ID:

Q96JM2

Full Name:

zinc finger protein 462

Calculated MW:

2506 aa, 285 kDa

Observed MW:

285 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB: 1:500-1:1000

IP: 0.5-4.0 µg for 1.0-3.0 mg of total protein lysate

Applications

Tested Applications:

WB, IP, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse retina tissue,

IP : mouse eye tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Mouse retina tissue lysates were subjected to SDS PAGE followed by western blot with 23793-1-AP (ZNF462 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.



IP result of anti-ZNF462 (IP:23793-1-AP, 4ug; Detection:23793-1-AP 1:500) with mouse eye tissue lysate 3600 ug.