

For Research Use Only

ZG16B Polyclonal antibody

Catalog Number: 23803-1-AP



Basic Information

Catalog Number:

23803-1-AP

Size:

650 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20740

GenBank Accession Number:

BC009722

GeneID (NCBI):

124220

UNIPROT ID:

Q96DAO

Full Name:

zymogen granule protein 16 homolog
B (rat)

Calculated MW:

208 aa, 23 kDa

Observed MW:

23 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : human urine sample,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human urine lysates were subjected to SDS PAGE followed by western blot with 23803-1-AP (ZG16B antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.