

For Research Use Only

CNGA2 Polyclonal antibody, PBS Only

Catalog Number: 23823-1-PBS



Basic Information

Catalog Number:

23823-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19899

GenBank Accession Number:

BC126302

GeneID (NCBI):

1260

UNIPROT ID:

Q16280

Full Name:

cyclic nucleotide gated channel alpha

2

Calculated MW:

664 aa, 76 kDa

Observed MW:

80 kDa

Purification Method:

Antigen Affinity purified

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Cyclic nucleotide-gated (CNG) channels are crucial for visual and olfactory transductions (PMID: 12432397). Native CNG channels are heterotetramers composed of homologous A and B subunits. The olfactory CNG channels are composed of two CNGA2 subunits, one CNGA4 and one CNGB1b subunit, each containing a cyclic nucleotide-binding domain (PMID: 22786723). CNGA2, also named as CNCA, CNCA1 and CNCG2, is the primary subunit of the olfactory CNG channel. Odorant signal transduction is probably mediated by a G-protein coupled cascade using cAMP as second messenger. Activation of olfactory channel by cyclic nucleotides leads to a depolarization of olfactory sensory neurons.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 23823-1-AP (CNGA2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 23823-1-PBS in a different storage buffer formulation.