

For Research Use Only

TRA2B Polyclonal antibody

Catalog Number: 23832-1-AP

1 Publications



Basic Information

Catalog Number:

23832-1-AP

Size:

230 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20683

GenBank Accession Number:

BC005898

GeneID (NCBI):

6434

UNIPROT ID:

P62995

Full Name:

transformer 2 beta homolog
(Drosophila)

Calculated MW:

34 kDa

Observed MW:

34 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Shu-Qi Li	32405436	J Adv Res	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

