

For Research Use Only

# ANXA2R Polyclonal antibody

Catalog Number: 23838-1-AP **1 Publications**



## Basic Information

<b>Catalog Number:</b> 23838-1-AP	<b>GenBank Accession Number:</b> BC067873	<b>Purification Method:</b> Antigen Affinity purified
<b>Size:</b> 550 µg/ml	<b>GeneID (NCBI):</b> 389289	<b>Recommended Dilutions:</b> IHC 1:20-1:200
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q3ZCQ2	
<b>Isotype:</b> IgG	<b>Full Name:</b> chromosome 5 open reading frame 39	
<b>Immunogen Catalog Number:</b> AG20807	<b>Calculated MW:</b> 193 aa, 22 kDa	

## Applications

<b>Tested Applications:</b> IHC, ELISA	<b>Positive Controls:</b> IHC : human prostate cancer tissue, human brain tissue
<b>Cited Applications:</b> WB	
<b>Species Specificity:</b> human	
<b>Cited Species:</b> human	
<b>Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0</b>	

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Ranran Sun	33008587	Biochem Biophys Res Commun	WB

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissue slide using 23838-1-AP (ANXA2R antibody at dilution of 1:50) (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissue slide using 23838-1-AP (ANXA2R antibody at dilution of 1:50) (under 40x lens).