

For Research Use Only

# ANXA2R Polyclonal antibody

Catalog Number: 23838-1-AP

1 Publications



## Basic Information

Catalog Number:

23838-1-AP

Size:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20807

GenBank Accession Number:

BC067873

GeneID (NCBI):

389289

UNIPROT ID:

Q3ZCQ2

Full Name:

chromosome 5 open reading frame 39

Calculated MW:

193 aa, 22 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

IHC 1:20-1:200

## Applications

Tested Applications:

IHC, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

**Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (\*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

Positive Controls:

IHC : human prostate cancer tissue, human brain tissue

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Ranran Sun	33008587	Biochem Biophys Res Commun	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissue slide using 23838-1-AP (ANXA2R antibody at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissue slide using 23838-1-AP (ANXA2R antibody at dilution of 1:50 (under 40x lens).