

For Research Use Only

RLN1 Polyclonal antibody

Catalog Number: 23852-1-AP



Basic Information

Catalog Number:

23852-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20880

GenBank Accession Number:

BC005956

GeneID (NCBI):

6013

UNIPROT ID:

P04808

Full Name:

relaxin 1

Calculated MW:

185 aa, 21 kDa

Purification Method:

Antigen Affinity purified

Applications

Tested Applications:

IF, IHC, ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

