

MYOZ3 Polyclonal antibody

Catalog Number: 23880-1-AP

Basic Information

Catalog Number:

23880-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20839

GenBank Accession Number:

BC074870

GeneID (NCBI):

91977

UNIPROT ID:

Q8TDC0

Full Name:

myozenin 3

Calculated MW:

251 aa, 27 kDa

Observed MW:

33 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB 1:500-1:2000

IHC 1:20-1:200

Applications

Tested Applications:

IHC, WB, ELISA

Species Specificity:

human, mouse, rat

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : rat skeletal muscle tissue, mouse skeletal muscle tissue

IHC : human skeletal muscle tissue,

Background Information

Storage

Storage:


Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

Selected Validation Data



rat skeletal muscle tissue were subjected to SDS PAGE followed by western blot with 23880-1-AP (MYOZ3 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human skeletal muscle tissue slide using 23880-1-AP (MYOZ3 Antibody) at dilution of 1:50 (under 40x lens).