

For Research Use Only

# RBM47 Polyclonal antibody

Catalog Number: 23902-1-AP

Featured Product

3 Publications



## Basic Information

Catalog Number:

23902-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20984

GenBank Accession Number:

BC126261

GeneID (NCBI):

54502

UNIPROT ID:

A0AV96

Full Name:

RNA binding motif protein 47

Calculated MW:

593 aa, 64 kDa

Observed MW:

64 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IP

Species Specificity:

human, rat, mouse

Cited Species:

human, mouse

Positive Controls:

WB : A549 cells, rat lung tissue

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Tao Guo	35831298	Cell Death Discov	WB
Pan He	31080349	Korean J Physiol Pharmacol	WB
Ziyao Li	38553613	Apoptosis	WB,IP

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



A549 cells were subjected to SDS PAGE followed by western blot with 23902-1-AP (RBM47 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.