

For Research Use Only

RBM47 Polyclonal antibody

Catalog Number: 23902-1-AP

Featured Product

2 Publications



Basic Information

Catalog Number:

23902-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20984

GenBank Accession Number:

BC126261

GeneID (NCBI):

54502

UNIPROT ID:

A0AV96

Full Name:

RNA binding motif protein 47

Calculated MW:

593 aa, 64 kDa

Observed MW:

64 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, rat, mouse

Cited Species:

human, mouse

Positive Controls:

WB : A549 cells, rat lung tissue

Background Information

Notable Publications

| Author | Pubmed ID | Journal | Application |
|---------|-----------|----------------------------|-------------|
| Tao Guo | 35831298 | Cell Death Discov | WB |
| Pan He | 31080349 | Korean J Physiol Pharmacol | WB |

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A549 cells were subjected to SDS PAGE followed by western blot with 23902-1-AP (RBM47 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.