

For Research Use Only

TMEM150 Polyclonal antibody

Catalog Number: 23912-1-AP



Basic Information

Catalog Number:

23912-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21019

GenBank Accession Number:

BC050466

GeneID (NCBI):

129303

UNIPROT ID:

Q86TG1

Full Name:

transmembrane protein 150

Calculated MW:

271 aa, 29 kDa

Observed MW:

33 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : Jurkat cells, mouse colon tissue, rat colon tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 23912-1-AP (TMEM150 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.