

For Research Use Only

ZBBX Polyclonal antibody

Catalog Number: 23958-1-AP



Basic Information

Catalog Number: 23958-1-AP	GenBank Accession Number: BC034229	Purification Method: Antigen Affinity purified
Source: Rabbit	GeneID (NCBI): 79740	Recommended Dilutions: WB: 1:500-1:3000 IHC: 1:50-1:500
Isotype: IgG	UNIPROT ID: A8MT70	
Immunogen Catalog Number: AG20988	Full Name: zinc finger, B-box domain containing	
	Calculated MW: 800 aa, 91 kDa	
	Observed MW: 88-96 kDa	

Applications

Tested Applications:

WB, IHC, ELISA

Species Specificity:

human, rat

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : A549 cells, rat brain tissue

IHC : mouse testis tissue,

Background Information

ZBBX (Zinc finger B-box domain-containing protein 1) contains a classic B-box-type zinc-finger (B-box1) motif that adopts a C6H2 zinc-binding fold and is therefore predicted to bind zinc ions and mediate protein-protein or protein-nucleic acid interactions. ZBBX was implicated previously by a de novo deletion in a male ADHD proband, and separately by a CNV gain disrupting ZBBX (PMID: 21832240).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 23958-1-AP (ZBBX antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded mouse testis tissue slide using 23958-1-AP (ZBBX antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).