

For Research Use Only

ZBBX Polyclonal antibody, PBS Only

Catalog Number: 23958-1-PBS



Basic Information

Catalog Number:

23958-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20988

GenBank Accession Number:

BC034229

GeneID (NCBI):

79740

UNIPROT ID:

A8MT70

Full Name:

zinc finger, B-box domain containing

Calculated MW:

800 aa, 91 kDa

Observed MW:

88-96 kDa

Purification Method:

Antigen Affinity purified

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, rat

Background Information

ZBBX (Zinc finger B-box domain-containing protein 1) contains a classic B-box-type zinc-finger (B-box1) motif that adopts a C6H2 zinc-binding fold and is therefore predicted to bind zinc ions and mediate protein-protein or protein-nucleic acid interactions. ZBBX was implicated previously by a de novo deletion in a male ADHD proband, and separately by a CNV gain disrupting ZBBX (PMID: 21832240).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 23958-1-AP (ZBBX antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 23958-1-PBS in a different storage buffer formulation.