

For Research Use Only

MPEG1 Polyclonal antibody

Catalog Number: 23960-1-AP



Basic Information

Catalog Number:

23960-1-AP

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21055

GenBank Accession Number:

BC104997

GeneID (NCBI):

219972

UNIPROT ID:

Q2M385

Full Name:

macrophage expressed 1

Calculated MW:

716 aa, 79 kDa

Observed MW:

70 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 23960-1-AP (MPEG1 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.