

For Research Use Only

# TTL2 Polyclonal antibody, PBS Only

Catalog Number: 23988-1-PBS



## Basic Information

Catalog Number:

23988-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21143

GenBank Accession Number:

BC047411

GeneID (NCBI):

83887

UNIPROT ID:

Q9BWW7

Full Name:

tubulin tyrosine ligase-like family,  
member 2

Calculated MW:

592 aa, 67 kDa

Observed MW:

67 kDa

Purification Method:

Antigen Affinity purified

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

## Background Information

TTL2 belongs to the Tubulin Tyrosine Ligase-Like (TTL) family, binds to tubulin and catalyzes the polyglutamylation modification of tubulin, participating in the organization and dynamic regulation of the microtubule cytoskeleton.

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 23988-1-AP (TLL2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 23988-1-PBS in a different storage buffer formulation.