

BAHD1 Polyclonal antibody

Catalog Number: 23991-1-AP

Basic Information

Catalog Number:

23991-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21164

GenBank Accession Number:

BC022782

GeneID (NCBI):

22893

UNIPROT ID:

Q8TBEO

Full Name:

bromo adjacent homology domain
containing 1

Calculated MW:

780 aa, 85 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

Gene repression in eukaryotes in response to developmental or environmental signals is a complex multistep process requiring the concerted action of many cellular factors, including DNA-binding transcription factors and cofactors. It also involves chromatin components and chromatin-binding/-modifying proteins, which are implicated in the establishment and maintenance of a condensed chromatin state, also referred to as facultative heterochromatin. BAHD1 is involved in heterochromatin formation and associates with the inactive X chromosome (Xi) in female cells. [PMID:19666599]

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

