

# CLRN2 Polyclonal antibody

Catalog Number: 23994-1-AP

1 Publications

## Basic Information

Catalog Number:

23994-1-AP

Size:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21182

GenBank Accession Number:

BC127863

GeneID (NCBI):

645104

UNIPROT ID:

A0PK11

Full Name:

clarin 2

Calculated MW:

232 aa, 25 kDa

Observed MW:

25 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB 1:200-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

mouse

Positive Controls:

WB : mouse eye tissue,

## Background Information

Clarín 2 belongs to the clarín family of genes. It's a small integral membrane glycoproteins with four transmembrane domains.

## Notable Publications

Author	Pubmed ID	Journal	Application
Lucy A Dunbar	31448880	EMBO Mol Med	WB

## Storage

Storage:


Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

## Selected Validation Data



mouse eye tissue were subjected to SDS PAGE followed by western blot with 23994-1-AP (CLRN2 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.