

For Research Use Only

C7orf63 Polyclonal antibody

Catalog Number: 24001-1-AP



Basic Information

Catalog Number:

24001-1-AP

Size:

260 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21205

GenBank Accession Number:

BC132771

GeneID (NCBI):

79846

UNIPROT ID:

A5D8W1

Full Name:

chromosome 7 open reading frame 63

Calculated MW:

941 aa, 106 kDa

Observed MW:

110 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse

Positive Controls:

WB : NIH/3T3 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with 24001-1-AP (C7orf63 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.