

For Research Use Only

# TMEM201 Polyclonal antibody

Catalog Number: 24092-1-AP **1 Publications**



## Basic Information

<b>Catalog Number:</b> 24092-1-AP	<b>GenBank Accession Number:</b> BC137291	<b>Purification Method:</b> Antigen Affinity purified
<b>Size:</b> 300 µg/ml	<b>GeneID (NCBI):</b> 199953	<b>Recommended Dilutions:</b> WB 1:200-1:1000
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q55NT2	
<b>Isotype:</b> IgG	<b>Full Name:</b> transmembrane protein 201	
<b>Immunogen Catalog Number:</b> AG21360	<b>Calculated MW:</b> 666 aa, 72 kDa	
	<b>Observed MW:</b> 43-45 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB : mouse brain tissue,
<b>Cited Applications:</b> WB, IHC	
<b>Species Specificity:</b> human, mouse	
<b>Cited Species:</b> human	

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Desheng Chen	37373430	Int J Mol Sci	WB, IHC


## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 24092-1-AP (TMEM201 Antibody) at dilution of 1:100 incubated at room temperature for 1.5 hours.