

For Research Use Only

ZMYM2 Polyclonal antibody

Catalog Number: 24112-1-AP **1 Publications**



Basic Information

Catalog Number:

24112-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21370

GenBank Accession Number:

BC036372

GeneID (NCBI):

7750

UNIPROT ID:

Q9UBW7

Full Name:

zinc finger, MYM-type 2

Calculated MW:

1377 aa, 155 kDa

Observed MW:

155-180 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB: 1:500-1:2000

IP: 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

Applications

Tested Applications:

WB, IP, ELISA

Species Specificity:

human

Positive Controls:

WB : Daudi cells, NCCIT cells, Raji cells

IP : Raji cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Ting Wei	40747902	J Med Chem	

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24112-1-AP (ZMYM2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 24112-1-AP (ZMYM2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



IP result of anti-ZMYM2 (IP:24112-1-AP, 4ug; Detection:24112-1-AP 1:600) with Raji cells lysate 1320 ug.