

For Research Use Only

ZMYM2 Polyclonal antibody, PBS Only

Catalog Number: 24112-1-PBS



Basic Information

Catalog Number:

24112-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21370

GenBank Accession Number:

BC036372

GeneID (NCBI):

7750

UNIPROT ID:

Q9UBW7

Full Name:

zinc finger, MYM-type 2

Calculated MW:

1377 aa, 155 kDa

Observed MW:

155-180 kDa

Purification Method:

Antigen Affinity purified

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24112-1-AP (ZMYM2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24112-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 24112-1-AP (ZMYM2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24112-1-PBS in a different storage buffer formulation.



IP result of anti-ZMYM2 (IP:24112-1-AP, 4ug; Detection:24112-1-AP 1:600) with Raji cells lysate 1320 ug. This data was developed using the same antibody clone with 24112-1-PBS in a different storage buffer formulation.