

For Research Use Only

# SCGF Polyclonal antibody

Catalog Number: 24138-1-AP



## Basic Information

Catalog Number:

24138-1-AP

Size:

340 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21365

GenBank Accession Number:

BC005810

GeneID (NCBI):

6320

UNIPROT ID:

Q9Y240

Full Name:

C-type lectin domain family 11,  
member A

Calculated MW:

36 kDa

Purification Method:

Antigen Affinity purified

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

## Background Information

CLEC11A, also named as CLECSF3, LSLCL, SCGF and p47, stimulates the proliferation and differentiation of hematopoietic precursor cells from various lineages, including erythrocytes, lymphocytes, granulocytes and macrophages. CLEC11A acts synergistically with other cytokines, including IL-3, GCSF, GM-CSF and FLT3 ligand. It suppresses SCF-stimulated erythrocyte proliferation.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

