

For Research Use Only

Multi tags Polyclonal antibody

Catalog Number: 24147-1-AP



Basic Information

Catalog Number:

24147-1-AP

Size:

800 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21457

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Transfected HEK-293 cells,

Background Information

This antibody recognize multi tags include: GST tag, 6*HIS tag, DDDDK tag, MYC tag and V5 tag.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with 24147-1-AP (Multi tag Antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.