

For Research Use Only

FOXB1 Polyclonal antibody

Catalog Number: 24285-1-AP



Basic Information

Catalog Number:

24285-1-AP

Size:

900 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18079

GenBank Accession Number:

BC113710

GeneID (NCBI):

27023

UNIPROT ID:

Q99853

Full Name:

forkhead box B1

Calculated MW:

325 aa, 35 kDa

Observed MW:

35 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : HepG2 cells, mouse liver tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 24285-1-AP (FOX B1 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.



mouse liver tissue were subjected to SDS PAGE followed by western blot with 24285-1-AP (FOX B1 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.