

For Research Use Only

# FAM71F2 Polyclonal antibody

Catalog Number: 24286-1-AP



## Basic Information

Catalog Number:

24286-1-AP

Size:

1500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16874

GenBank Accession Number:

BC105730

GeneID (NCBI):

346653

UNIPROT ID:

Q6NXP2

Full Name:

family with sequence similarity 71, member F2

Calculated MW:

309 aa, 35 kDa

Observed MW:

40-45 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

IF/ICC 1:20-1:200

## Applications

Tested Applications:

WB, IF/ICC, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : PC-3 cells, mouse liver tissue

IF/ICC : PC-3 cells,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



PC-3 cells were subjected to SDS PAGE followed by western blot with 24286-1-AP (FAM71F2 Antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (-20°C Ethanol) fixed PC-3 cells using 24286-1-AP (FAM71F2 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).