

For Research Use Only

C7orf31 Polyclonal antibody

Catalog Number: 24354-1-AP



Basic Information

Catalog Number:

24354-1-AP

Size:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19530

GenBank Accession Number:

BC043269

GeneID (NCBI):

136895

UNIPROT ID:

Q8N865

Full Name:

chromosome 7 open reading frame 31

Calculated MW:

590 aa, 68 kDa

Observed MW:

65 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, Mouse, Rat

Positive Controls:

WB : mouse testis tissue, rat testis tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 24354-1-AP (C7orf31 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.