

For Research Use Only

GPATCH2 Polyclonal antibody, PBS Only

Catalog Number: 24366-1-PBS



Basic Information

Catalog Number:

24366-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21647

GenBank Accession Number:

BC063474

GeneID (NCBI):

55105

UNIPROT ID:

Q9NW75

Full Name:

G patch domain containing 2

Calculated MW:

528 aa, 59 kDa

Observed MW:

65-70 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



mouse spleen tissue were subjected to SDS PAGE followed by western blot with 24366-1-AP (GPATCH2 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24366-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human breast cancer tissue slide using 24366-1-AP (GPATCH2 Antibody) at dilution of 1:50 (under 40x lens). This data was developed using the same antibody clone with 24366-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human testis tissue slide using 24366-1-AP (GPATCH2 Antibody) at dilution of 1:50 (under 40x lens). This data was developed using the same antibody clone with 24366-1-PBS in a different storage buffer formulation.