

For Research Use Only

SETD6 Polyclonal antibody

Catalog Number: 24377-1-AP



Basic Information

Catalog Number: 24377-1-AP	GenBank Accession Number: BC022451	Purification Method: Antigen affinity purification
Size: 1100 µg/ml	GeneID (NCBI): 79918	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q8TBK2	
Isotype: IgG	Full Name: SET domain containing 6	
Immunogen Catalog Number: AG19485	Calculated MW: 473 aa, 53 kDa	
	Observed MW: 53 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : RAW 264.7 cells,
Species Specificity: human, mouse	

Background Information

SETD6 also named as SET domain containing 6 is a 473 amino-acid protein, which localizes in the Nucleus. The molecular weight of SETD6 is 53kDa. SETD6 consists of 1 SET domain and belongs to the class V-like SAM-binding methyltransferase superfamily. SETD6 as a protein-lysine N-methyltransferase that specifically monomethylates at Lysine 310 of the RELA subunit of NF-kappa-B complex, leading to down-regulate NF-kappa-B transcription factor activity.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with 24377-1-AP (SETD6 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.