

SETD6 Polyclonal antibody

Catalog Number: 24377-1-AP

Basic Information

Catalog Number:

24377-1-AP

Size:

1100 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19485

GenBank Accession Number:

BC022451

GeneID (NCBI):

79918

UNIPROT ID:

Q8TBK2

Full Name:

SET domain containing 6

Calculated MW:

473 aa, 53 kDa

Observed MW:

53 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : RAW 264.7 cells,

Background Information

SETD6 also named as SET domain containing 6 is a 473 amino-acid protein, which localizes in the Nucleus. The molecular weight of SETD6 is 53kDa. SETD6 consists of 1 SET domain and belongs to the class V-like SAM-binding methyltransferase superfamily. SETD6 as a protein-lysine N-methyltransferase that specifically monomethylates at Lysine 310 of the RELA subunit of NF-kappa-B complex, leading to down-regulate NF-kappa-B transcription factor activity.

Storage

Storage:


Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

Selected Validation Data



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with 24377-1-AP (SETD6 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.