

For Research Use Only

METTL13 Polyclonal antibody, PBS Only

Catalog Number: 24379-1-PBS



Basic Information

Catalog Number: 24379-1-PBS	GenBank Accession Number: BC029083	Purification Method: Antigen affinity purification
Concentration: 1 mg/ml	GeneID (NCBI): 51603	
Source: Rabbit	UNIPROT ID: Q8N6R0	
Isotype: IgG	Full Name: methyltransferase like 13	
Immunogen Catalog Number: AG19355	Calculated MW: 699 aa, 79 kDa	
	Observed MW: 72-78 kDa	

Applications

Tested Applications:
WB, IF/ICC, Indirect ELISA

Species Specificity:
human

Background Information

Protein methyltransferase like 13 (METTL13), also called FEAT, is coded by gene METTL13, which is located at chromosome 1q24.3. Purified from rat livers in 2011 by a Japanese researcher, METTL13 was found to be abnormally overexpressed in several human cancers including lung cancer and to drive tumorigenesis in transgenic mice. A study indicated that METTL13 inhibits apoptosis of lung, breast and liver cancer cells and miR-16 can repress the expression of METTL13 by binding to its 3' -untranslated region. (PMID: 33985542)

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24379-1-AP (METTL13 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24379-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed U-251 cells using METTL13 antibody (24379-1-AP) at dilution of 1:400 and CoraLite® 488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 24379-1-PBS in a different storage buffer formulation.