

For Research Use Only

MPND Polyclonal antibody

Catalog Number: 24420-1-AP



Basic Information

Catalog Number:

24420-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19531

GenBank Accession Number:

BC032652

GeneID (NCBI):

84954

UNIPROT ID:

Q8N594

Full Name:

MPN domain containing

Calculated MW:

471 aa, 51 kDa

Observed MW:

49-55 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : SMMC-7721 cells, HepG2 cells, mouse liver tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



SMMC-7721 cells were subjected to SDS PAGE followed by western blot with 24420-1-AP (MPND Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.